

FCP-06 AUTO CHART VISION

There is an SD card behind the monitor, update the software by it.

Software and information are retained when power is turned off.

Four units can be operated at the same time without problem.

Four kinds of astigmatism clock dial:  
①broken kine 15° ②real line 15°  
③broken line 10° ④real line 10°

Operation and controls are simple and easy by remote contorl , can link with Auto phoropter.

No copyright issue.

Software are all developed by ourselves.

Specification of screen

- 1. Panel:19"(4:3)
- 2. Pixel pitch: 0.285mm
- 3. Max.brightness:300cd/m2
- 4. Max.contrast ratio: 800:1
- 5. Max.resolution:1280×1024
- 6. Input signal: RGB
- 7. Sigal joint: 15Pin D Pattern Plug
- 8. Energy consumption: 35W
- 9. Net weight: 5KG
- 10. Compatibility: PC/MAC
- 11. Dimension:
- 60 (D) ×425(W) ×440(H)mm

Placement

- 1. Can be hung on the wall by a bracket.
- 2. Can stand on ground by a support.
- 3. Can be set on the table.

FCP-06

Crescent Health Co., Ltd


ADD: Rm 14-1002,No.388, Wanrong Rd, Zhabei,Shanghai,China

TEL:0086-021-62546941 FAX:0086-021-32585057

http://www.crescenthealthco.com









Email:Brook822@163.com & Brook@crescenthealthco.com

COLOFUL SCREEN





Offers 19" panel with high clarity chart projection.  
And accurate examination,  
Moves with no noise.

FCP-06 provides single letter and horizontal /vertical line masks, it can meet all the function examination.



SELECTABLE WORKING DISTANCE



working distance : (2-6.5m)  
The chart size automatically  
calibrates to the working distance.  
the suitable working distance: (3-6m)

TWO TYPE:

Chart Type A




Chart Type B

