


Crescent Health Co., Ltd

Ophthalmic Surgical Instruments



Crescent Health Co., Ltd

ADD: Rm 14-1002, No. 388, Wanrong Rd,

Zhabei, Shanghai, China

TEL: 0086-021-62546941 FAX: 0086-21-52662361

EMAIL: info@crescenthealthco.com



www.crescenthealthco.com

CONTENT

新品	New Surgical Instruments	1-18
20G新型眼内器械	20G New Inteaoocular Medical Instruments	1-9
新型斜视手术器械	New Steabismus Surgical Instruments	10-12
新材料手术器械	New Material Surgical Instruments	13-14
镀金刚石膜手术器械	Resembling Diamond Film Layer	14-16
镀金手术器械	New Surgical Instruments With The Technology Of Titanium Nitride	17-18
宝石手术刀	Sapphire Surgical Knife	18-19
黑钻石手术刀	Black diamond Surgical Knife	18
蓝宝石手术刀	Sapphire Surgical Knife	19
开睑器	Eye Speculum	20-22
钩	Hook	23-24
移核器械	Manipulation Instruments	24
刀	Knife	25-27
切口刀	Cut Knife	25-26
劈核刀	Cystotome Knife	27
复位分离器	Spatula	28-29
匙	Spoon	29
白内障匙	Lens Spoop	29
睑板腺囊肿锐匙	Chalazion Spoon	29
晶体线环	Lens Loop	30
普通线环	Ordinary Loop	30
注水式线环	Infusion Lens Loop	30
镊子	Forceps	31-43

晶体植入镊/植入器	IOL Implantation Forceps / Holder	31-32
撕囊镊	Capsulorhexis Forceps	32-34
缝线结扎镊	Suturing Forceps	34-36
系结镊	Tying Forceps	36-39
无损伤镊	Corneal Scleral Forceps	39
其他普通镊类（虹膜镊、眼用镊、结合膜镊、异物镊、固定镊、斜视镊、囊肿镊、沙眼镊、翻眼镊、睫毛镊）	Others (Iris Forceps, Ophthalmic Forceps, Foreign Body Forceps, Fixation Forceps)	40-43

剪刀 Scissors 43-47

角膜剪、小梁剪	Corneal Scissors, Micro-corneal Scissors	43-44
膜状内障剪	Cataract Scissors	44-45
维纳斯剪	Vannas Scissors	45-46
眼用手术剪	Ophthalmic Scissors	46
其他普通剪刀（虹膜剪、切腱剪、巩膜剪、斜视剪、眼球摘出剪）	Others (Iris Scissors, Sclera Scissors, Strabismus Scissors)	47-48

持针器 Holders 47-48

注吸针 Needles 49-52

注吸器	Cannula	50-51
冲洗针	Irrigation Needle	51

屈光手术器械 Diopter Instruments 52-53

其他手术器械 Others Instruments 54-60

手术器械包 Surgical Instrument Packages 61-71

The meaning of letter of the products name

A: Improved; B: Improved; C: Improved; D: Special Material

T: Titanium; Z: Black Diamond Sapphire Surgical Knife

S: Sapphire Surgical Knife; M: Plating Diamond Film Layer

X: The Head can be rotated;

H: Not only use new material and plating diamond technology, but also apply the latest technology of titanium nitride in the handle which makes it to be yellow.

20G新型眼内器械

20G New Intraocular Instruments



产品采用进口高强度不锈钢材料，具有硬度高、耐腐蚀和抗酸性的特点；镀金刚石膜工艺，增强了刃口锋利度，提高耐磨耐腐蚀性，减少眩光、降低视觉疲劳；特种焊接工艺，即使是高温高压消毒后依然伸缩自如；其中旋转式眼内器械头部360°旋转，满足各种角度手术要求；彩色标识，术中快速分辨器械类别，提高手术效率；头部可拆卸，配置专用清洗工具，清洗更方便。

The product is made of imported alloy steel with superior intension, which is stainless and nondeformable. Resembling diamond film layer enhances the sharpness of cutting edge, improves the abrasionproof and corrosion proof and reduces the glare and vision tiredness. The products still moves freely after sterilization under high temperature and high pressure, which benefits from the special welding technology.

旋转式眼内镊 Rotation Vireoretinal Forceps



53190TX
等角镊
Equiangular Forceps

等角型，头部0.25×0.25mm
Equiangular, Tip 0.25×0.25mm



53191TX
末端夹持镊
Retinal Forceps

微型头设计，头部0.15×0.4mm
Micro tip, Tip 0.15×0.4mm



53568TX
内界膜镊
Membrana Limitans
Internia Forcep



直型、头部0.3×0.2mm
Straight, Tip 0.3×0.2mm



53185TX
弯形镊
Curved Forcep



头部0.2×0.15mm
Tip 0.2×0.15mm



53590TX
直镊
Straight Forcep



尖头形、头部0.15×0.2mm
Tip 0.15×0.2mm



53178TX
钉形镊
Spiked Forcep



头部0.15×0.15mm
Tip 0.15×0.15mm



53192TX
网膜下镊
Retinal Forcep



头部0.2×0.25mm
Tip 0.2×0.25mm



53174TX
三爪异物镊
Three-claws Foreign Body Forcep



三爪, ϕ 0.2mm
Three claws, ϕ 0.2mm



53184TX
杯形夹口镊
Double Cup Forcep



头部1×0.4mm
Tip 1×0.4mm

旋转式眼内剪 Rotation Vireoretinal Scissors



54290TXM
水平剪
Horizontal Scissor



镀类金刚石膜, 矛型, 刃长2.4mm
Resembling diamond film layer, spear shaped, 2.4mm in blade length.



54230TXM
垂直剪
Vertical Scissor


镀类金刚石膜, 钢刀直型, 刃长1.55mm
Resembling diamond film layer, 1.55mm in blade length.



54289TXM
弯形剪
Curved Scissor


镀类金刚石膜, 月弯型, 刃长1.3mm
Resembling diamond film layer, falcated, 1.3 mm in blade length.

三片式眼内镊 Three Piece Type of Vireoretinal Forceps




53172TS
三片式异物镊
Ophthalmic Foreign Body Forceps

无齿45° 全长148mm
No teeth 45°, 148mm in total length



53173TS
三片式异物镊
Ophthalmic Foreign Body Forceps

网纹齿45° 全长148mm
Textured teeth 45°, 148mm in total length



53174TS
三片式异物镊
Ophthalmic Foreign Body
Forceps (for Vitreous Body)

三爪0.2×0.5 mm 全长148mm
Three claws, 148mm in total length



53175TS
三片式视网膜镊
Micro-retina Forceps No teeth

无齿1×0.4 mm 全长148mm
1×0.4 mm, 148mm in total length



53176TS
三片式视网膜镊
Micro-retina Forceps Straight list teeth

直纹齿1×0.4 mm 全长148mm
1×0.4 mm, 148mm in total length



53178TS
三片式眼内镊
Intraocular Forcep

钉形 0.2×0.25mm 全长148mm
Spiked, 0.2×0.25mm 148mm in total length





53185TS
三片式眼内镊
Intraocular Forcep



弯圆角0.2×0.15mm全长148mm
Tip 0.2×0.15mm ,148mm in total length



53568TS
三片式内界膜镊
Membrana Limitans Interna Forcep Straight,



直型0.3×0.2mm 全长148mm
0.3mm in tip width,148mm in total length



53590TS
三片式眼内镊
Intraocular Forcep



尖头型0.2×0.25mm 全长148mm
Tip 0.2×0.25mm ,148mm in total length

三片式眼内剪 Three Piece Type of Vireoretinal Scissor



54289TSM
三片式眼内剪
Intraocular Scissor



月弯型刃2.8mm刃长, 全长148mm
Falcated ,2.4mm in blade length ,148mm in total length



54290TSM
三片式眼内剪
Intraocular Scissor



矛型刃 2.4mm刃长, 全长148mm
Spear Shaped ,2.4mm in blade length ,148mm in total length

23G眼内器械 23G Intraocular Instruments



53140TS
三片式眼内镊
Intraocular Forcep



23G 尖头型0.15×0.2mm 全长148mm
23G Shaped Tip 0.2×0.2mm ,148mm in total length



53148TS
三片式眼内镊
Intraocular Forcep



23G 钉形 0.15×0.15mm 全长148mm
23G Spiked ,0.15×0.15mm 148mm in total length



53569TS
三片式内界膜镊



23G 直型 内界膜镊0.25×0.15 网状 全长148mm



53186TS
三片式眼内镊
Intraocular Forcep

23G 弯形镊 0.8×0.15mm 网状 全长148mm



54291TSM
三片式眼内剪
Intraocular Scissor

23G 月弯形R6 刃长1.6mm 全长148mm

飞秒激光手术器械

Femtosecond laser surgery Instruments

国内首套专门为飞秒激光手术精心设计的手术器械，产品采用新材料新工艺精制而成，品质卓越，是飞秒激光手术的最佳选择。

The product is firstly and specially designed by both 66 Vision Tech Co.,Ltd and the national top specialist in china, which is made of new material with the new technology and have good quality. It is the best choice for Femtosecond laser surgery.

虹膜复位器 Iris Spatula



52429T

双头式，一头椭圆0.2×0.4mm，一头扁圆1.2×0.5mm 全长140mm，
Double-headed, Head width of 0.2×0.4mm and 1.2×0.5mm, 140mm in total length

飞秒透镜分离器 Femtosecond laser surgery lens separator



52395T

双头，左端φ0.3mm，右端φ0.2mm
Double-headed, φ0.3mm and φ0.2mm

飞秒分离匙 Femtosecond laser surgery corneal flap separator



52435T

头端φ0.3mm
φ0.3mm tip

飞秒透镜镊 Femtosecond laser surgery lens Forcep



53155T

弯形（网纹）0.2×0.15mm
Curved Textured 0.2×0.15mm

系结镊 Tying Forceps



53338DH

直形、全长110mm
Straight, 110mm in total length

新型斜视手术器械

New Strabismus Surgical Instruments

国内首套专门为斜视手术精心设计的斜视手术器械；与国内首席斜弱视专家联合开发；产品采用新材料新工艺加工而成，品种齐全，品质卓越，是斜视手术最完美的选择！

The product is firstly and specially designed by both 66 Vision Tech Co., Ltd and the national top strabismus specialist for the strabismus surgery in China, which is made of new material with the new technology and has various models and good quality. It is your best choice for the strabismus surgery.

斜视镊 Micro-corneal Scissors



53800A

左式弯平台、0.25唇头钩，全长105mm
Left curved platform, 1×2 teeth of 0.25mm 105mm in total length

53801A

右式弯平台、0.25唇头钩，全长105mm
Right curved platform, 1×2 teeth of 0.25mm 105mm in total length

53802A

直式平台、0.25唇头钩，全长105mm
Straight platform, 1×2 teeth of 0.25mm 105mm in total length

斜视镊 Micro-corneal Scissors



53583

直式平台、0.35唇头钩，全长88mm
Straight platform, 1×2 teeth of 0.35mm 88mm in total length

斜视钩 Strabismus Hook



52456A

钝钩头，头宽2.5mm，全长121mm
Dull tip, 2.5mm 121mm in total length

斜视钩 Strabismus Hook



52459A

F钩头，头宽1.5mm，全长129mm
F head, 1.5mm 129mm in total length

斜视钩 Strabismus Hook



52439A

弯圆头，头宽1.8mm，全长130mm
Curved round tip, 1.8mm in tip width 130mm in total length

52440

弯圆头，头宽1.6mm，全长129mm
Curved round tip, 1.6mm in tip width 129mm in total length



52457A

立体钩，椭圆头0.7mm，全长125mm
Tridimensional hook, elliptical head 0.7mm, 125mm in total length

斜视钩 Strabismus Hook



52451A

弯圆头，头宽1.2mm，全长126mm
Curved round tip, 1.2mm in tip width 126mm in total length

52452A

弯圆头，头宽1.3mm，全长126mm
Curved round tip, 1.3mm in tip width 126mm in total length

52455A

弯圆头，头宽0.6mm，全长121mm
Curved round tip, 0.6mm in tip width 121mm in total length

斜视钩 Strabismus Hook



52458A

刻度头，头宽1.75mm，全长116mm
Head with dial, 1.75mm in tip width 116mm in total length

开睑器 Eye Speculum



51280

遮瞳式，板宽6mm，扩张杆长11.5mm，全长64mm
Solid Blades, 6mm in board width, 11.5mm blades in length, 64mm in total length

51280A

遮瞳式，板宽8mm，扩张杆长14mm，全长64mm
Solid Blades, 8mm in board width, 14mm blades in length, 64mm in total length

斜视钩 Strabismus Hook



52460A

双钩头、直径0.8mm, 全长119mm
SDouble head, $\varnothing$ 0.8mm 119mm in total length

斜视圆挠板 Strabismus Hook



52462A

头宽9.2, 全长118mm
9.2mm in tip width 118mm in total length

52463

头宽6mm, 全长127mm
6mm in tip width 127mm in total length

眼睑拉钩 Eyelid Hook



51452A

头宽10mm, 全长121mm
10mm in tip width 121mm in total length

显微持针钳 Micro Needle Holder



55159T

弯锁止钩, 全长135mm
Curved with lock, 135mm in total length

新材料手术器械

New Material

新材料手术器械采用进口高强度不锈钢新材料精制而成, 产品具有硬度高、耐腐蚀和抗酸性的特点。

New material surgical instruments adopt imported high strength stainless steel refined but become, product has high hardness, corrosion resistance and acid resistance, etc

撕囊镊 Capsulorhexis Forceps



53402D

弧形角形头, 全长110mm
Curved, angled, 110mm in total length



53419D

弧形瓜子头, 全长110mm
Curved, 1-claw-shaped head, 110mm in total length

缝线结扎镊 Suturing Forceps



53411D

直平台, 唇头钩0.15, 全长110mm
Straight, platform, 1 $\times$ 2 teeth of 0.15mm, 110mm in total length

53412D

直平台, 唇头钩0.20, 全长110mm
Straight, platform, 1 $\times$ 2 teeth of 0.20mm, 110mm in total length

53413D

直平台, 唇头钩0.12, 全长110mm
Straight, platform, 1 $\times$ 2 teeth of 0.12mm, 110mm in total length

人工晶体植入镊 IOL Implantation Forceps



53401D

角形, 头部直径 $\varnothing$ 0.6mm, 全长110mm
Angled, $\varnothing$ 0.6mm of tip, 110mm in total length

系结镊 Tying Forceps



53414D

直平台, 全长110mm
Straight platform, 110mm in total length



53415D

弯平台, 全长110mm
curved platform,



53416D

长角形平台, 全长110mm
Long angled platform,

53417D

短角形平台, 全长110mm
Short angled platform



角膜上皮镊 Corneal Epidermic Forceps



53418D

圆形, $\phi 0.6\text{mm}$, 全长110mm
Round, $\phi 0.6\text{mm}$, 110mm in total length



缝线结扎镊 Suturing Forceps



53670D

大齿, 0.12mm, 全长109mm
0.12mm, 109mm in total length

53671D

大齿, 0.3mm, 全长109mm
0.3mm, 109mm in total length



镀金刚石膜手术器械

Resembling Diamond Film Layer

镀金刚石膜是指采用特殊工艺在手术器械头部表面镀一层薄膜。膜层具有高硬度、耐磨损和耐腐蚀的特殊性能, 锋利度明显提高, **剪切寿命提高5-10倍以上**。深褐色膜层能防止眩光, **有效降低视觉疲劳**。膜层极其光滑, 具有良好生物相容性, 并通过人性化设计, 医生可安全而极其顺畅地开展高质量的手术。

Resembling diamond film adopts special technology which the head of surgical instrument covering with a film. The film has the superior performances of high-rigidity, abrasionproof and corrosion proof. The sharpness of scissors and cutting life of medical instrument which **covering film will increased over 5 to 10 times**. The puce film can prevent the glare and **reduce the vision tiredness effectively**. In addition, the smoothness of film and biological compatibility, doctors can operate safely and smoothly.

钛合金镀金刚石膜手术器械

Titanium Material with Resembling Diamond Film Layer

角膜剪 Corneal Scissors



54109TM

镀类金刚石膜, 弯形, 刃长14mm, 全长120mm
Plating diamond film layer, slight curved, 14mm in blade length, 120mm in total length



54053TM

镀类金刚石膜, 直形, 刃长14mm, 全长120mm
Plating diamond film layer, straight, 14mm in blade length, 120mm in total length

小梁剪 Micro-corneal Scissors



54108TM

镀类金刚石膜, 弯尖头, 刃长12mm, 全长120mm
Plating diamond film layer, curved tips, 12mm in blade length, 120mm in total length



54110TM

镀类金刚石膜, 直尖头, 刃长12mm, 全长120mm
Plating diamond film layer, straight tips, 12mm in blade length, 120mm in total length

膜状内障剪 Cataract Scissors Membrane



54268TM

镀类金刚石膜, 直形, 刃长16mm, 全长122mm
Plating diamond film layer, straight, 16mm in blade length, 122mm in total length



54269TM

镀类金刚石膜, 弯形, 刃长16mm, 全长122mm
Plating diamond film layer, slight curved, 16mm in blade length, 122mm in total length



新材料镀金刚石膜手术器械

New Material with Resembling Diamond Film Layer

角膜剪 Corneal Scissors



54109BDM

镀类金刚石膜，弯形，刃长14mm，全长120mm
Plating diamond film layer, slight curved, 14mm in blade length, 120mm in total length

小梁剪 Micro-corneal Scissors



54108BDM

镀类金刚石膜，弯尖，刃长12mm，全长120mm
Plating diamond film layer, curved, 12mm in blade length, 120mm in total length



54110BDM

镀类金刚石膜，直尖，刃长12mm，全长120mm
Plating diamond film layer, straight tips, 12mm in blade length, 120mm in total length

膜状内障剪 Cataract Scissors Membrane



54268BDM

镀类金刚石膜，直形，刃长16mm，全长120mm
Plating diamond film layer, straight, 16mm in blade length, 120mm in total length



54269BDM

镀类金刚石膜，弯形，刃长16mm，全长120mm
Plating diamond film layer, slight curved, 16mm in blade length, 120mm in total length

镀金手术器械

New Surgical Instruments With The Technology Of Titanium Nitride

镀金手术器械不仅采用了新材料新工艺镀金膜技术，而且在手柄部采用了最新的氮化钛涂层工艺，使手柄呈金黄色，增强了产品的视觉效果和操作舒适感。

The surgical instruments not only use the new material and resembling diamond technology, but also apply the latest technology of titanium nitride in the handle which makes it to be yellow.

缝线结扎镊 Suturing Forceps



53411DH

直平台，唇头钩0.15，全长110mm
Straight, platform, 1×2 teeth of 0.15mm, 110mm in total length

53412DH

直平台，唇头钩0.20，全长110mm
Straight, platform, 1×2 teeth of 0.20mm, 110mm in total length

53413DH

直平台，唇头钩0.12，全长110mm
Straight, platform, 1×2 teeth of 0.12mm, 110mm in total length

系结镊 Tying Forceps



53414DH

直平台，全长110mm
Straight platform, 110mm in total length



53415DH

弯平台，全长110mm
curved platform, 110mm in total length



53416DH

长角形平台，全长110mm
Long angled platform, 110mm in total length

53417DH

短角形平台，全长110mm
Short angled platform, 110mm in total length

角膜剪 Corneal Scissors



54109BHM

镀类金刚石膜，弯形，刃长14mm，全长120mm
Plating diamond film layer, slight curved, 14mm in blade length, 120mm in total length

54053BHM

镀类金刚石膜，直形，刃长14mm，全长120mm
Plating diamond film layer, straight, 14mm in blade length, 120mm in total length

小梁剪 Micro-corneal Scissors



54108BHM

镀类金刚石膜，弯尖，刃长12mm，全长120mm
Plating diamond film layer, curved, 12mm in blade length, 120mm in total length

新型宝石手术刀 Surgical Knife

刀片：黑钻石 Blade: Black diamond
刀柄：钛合金 Handle: Titanium alloy



50102ZA

月形隧道刀 2.8×0.25mm
Tunnel knife 2.8×0.25mm



50401ZA

弯头穿刺刀 2.8×0.2mm×60°
Curved, 2 cutting edges 2.8×0.2mm×60°



50403ZA

弯头穿刺刀 3.0×0.2mm×60°
Curved, 2 cutting edges 3.0×0.2mm×60°



50405ZA

弯头穿刺刀 3.2×0.2mm×60°
Curved, 2 cutting edges 3.2×0.2mm×60°



50406ZA

弯头穿刺刀 2.65×0.2mm×60°
Curved, 2 cutting edges 2.65×0.2mm×60°



50409ZA

弯头穿刺刀 2.2mm×0.2mm×60°
Curved, 2 cutting edges 2.2mm×0.2mm×60°



50502ZA

直头穿刺刀 1.0×0.25mm×60°
Straight, 2 cutting edges 1.0×0.25mm×60°



50500ZA

直头穿刺刀 1.0×0.25mm×35°
Straight, 2 cutting edges 1.0×0.25mm×35°



50200ZA

双面刃切口刀 1.0×0.2mm×30°
Bi-facet blade 1.0×0.2mm×30°



50205ZA

双面刃切口刀 1.0×0.2mm×45°
Bi-facet blade 1.0×0.2mm×45°

新型宝石手术刀 Sapphire Surgical Knife

刀片：蓝宝石 Blade: Sapphire
刀柄：钛合金 Handle: Titanium alloy



50102TA

月形隧道刀 2.8×0.25mm
Tunnel knife 2.8×0.25mm



50401TA

弯头穿刺刀 3.0×0.2mm×60°
Curved, 2 cutting edges 3.0×0.2mm×60°



50403TA

弯头穿刺刀 3.0×0.2mm×60°
Curved, 2 cutting edges 3.0×0.2mm×60°



50405TA

弯头穿刺刀 3.2×0.2mm×60°
Curved, 2 cutting edges 3.2×0.2mm×60°



50200TA

双面刃切口刀 1.0×0.2mm×30°
Bi-facet blade 1.0×0.2mm×30°



50205TA

双面刃切口刀 1.0×0.2mm×45°
Bi-facet blade 1.0×0.2mm×45°



50500TA

直头穿刺刀 1.0×0.25mm×35°
Straight, 2 cutting edges 1.0×0.25mm×35°



50502TA

直头穿刺刀 1.0×0.25mm×60°
Straight, 2 cutting edges 1.0×0.25mm×60°



50600TA

单面刃切口刀 1.0×0.2mm×30°
Single-facet blade 1.0×0.2mm×30°



50901TA

三角四刃刀, 3.5×5.16mm
Four Blade 3.5×5.16mm

开睑器 Eye Speculum



51262

露睫式、扩张杆长11mm、全长75mm
Open blades, stainless 11mm blade in length,
75mm in total length

开睑器 Eye Speculum



51272A

遮睫式、扩张杆长15mm、全长75mm
Solid blades, stainless 15mm blade in length,
75mm in total length

开睑器 Eye Speculum



51281A

露睫式、扩张杆长15mm、全长44mm
Open blades, stainless 15mm blade in length,
44mm in total length

开睑器 Eye Speculum



51282A

遮睫式、扩张杆长14mm、全长46mm
Solid blades, stainless 14mm blade in length,
46mm in total length

开睑器 Eye Speculum



51283A

露睫式、扩张杆长14mm、全长64mm
Open blades, stainless 14mm blade in length,
64mm in total length

51283T

露睫式、扩张杆长14mm、全长64mm
Open blades, stainless 14mm blade in length,
64mm in total length

开睑器 Eye Speculum



51284T

露睫式、扩张杆长12mm、全长70mm
Open blades, stainless 12mm blade in length,
70mm in total length

开睑器 Eye Speculum



51285T

露睫式、扩张杆长12mm、全长70mm
Open blades, stainless 12mm blade in length,
70mm in total length

开睑器 Eye Speculum



51286A

钢丝式、扩张杆长14mm、全长50mm
Wire speculum, stainless 14mm blade in length,
50mm in total length

开睑器 Eye Speculum



51288A

钢丝式、扩张杆长15mm、全长50mm
Wire speculum, stainless 15mm blade in length,
50mm in total length


开睑器 Eye Speculum



51290A

钢丝式、扩张杆长11mm、全长37mm
Wire speculum, stainless 11mm blade in length,
37mm in total length

开睑器 Eye Speculum



51295

扩张杆长14mm、全长65mm
Stainless 14mm blade in length, 65mm in total length

开睑器 Eye Speculum



51302

速捷式、扩张杆长14mm、全长90mm
Solid blades, Stainless 14mm blade in length, 90mm in total length

眼睑拉钩 Eyelid Retractor



51450

头宽8mm、全长125mm
8mm in blade width, 125mm in total length



51452

头宽10mm、全长130mm
10mm in blade width, 130mm in total length

51454

头宽12mm、全长135mm
12mm in blade width, 135mm in total length

人工晶体定位钩 IOL Orientation Hook



51474A

角形 $\phi 0.2\text{mm}$ 、全长110mm
Angled $\phi 0.2\text{mm}$, 110mm in total length

51474T

角形 $\phi 0.2\text{mm}$ 、全长110mm
Angled $\phi 0.2\text{mm}$, 110mm in total length

人工晶体定位钩 IOL Orientation Hook



51478A

圆弯 $\phi 0.2\text{mm}$ 、15° 全长110mm
0.2mm cylindrical tip, angled 15°,
110mm in total length

51478T

圆弯 $\phi 0.2\text{mm}$ 、15° 全长110mm
0.2mm cylindrical tip, angled 15°,
110mm in total length



51479T

圆弯 $\phi 0.2\text{mm}$ 、40° 全长110mm
0.2mm cylindrical tip, angled 40°,
110mm in total length

斜视钩 Strabismus Hook



52438A

带孔 $\Phi 0.6$ 全长130mm
Head aperture $\Phi 0.6$ mm, 130mm in total length



斜视钩 Strabismus Hook



52450

弯圆头、头宽8mm、全长120mm
Blunt, curved tip, 8mm in head width, 120mm in total length

52452

弯圆头、头宽10mm、全长130mm
Blunt, curved tip, 10mm in head width, 130mm in total length

52454

弯圆头、头宽12mm、全长135mm
Blunt, curved tip, 12mm in head width, 135mm in total length



核转动器 Lens Manipulator



51470A

叉形长头、头宽0.7mm、全长115mm
"Y" shaped tip, 0.7mm in head width, 115mm in total length

51470T

叉形长头、头宽0.7mm、全长115mm
"Y" shaped tip, 0.7mm in head width, 115mm in total length



移核器 Lens Nucleus Hook



51471A

弯角形、头宽0.2mm、全长115mm
Curved and angled, 0.2mm in head width, 115mm in total length



前房角切开刀 Goniotomy Knife



52231A

直形、刃长2mm、全长115mm
Straight, 2mm in blade length, 115mm in total length



隧道刀 Discission Knife



52274A

尖头、宽2.1mm、全长128mm
Tine, 2.1mm in width, 128mm in total length



隧道刀 Discission Knife



52275A

圆头、宽2.5mm、全长128mm
Round Tip, 2.5mm in width, 128mm in total length



碎核刀 Cataract Knife



52276D

角形、刃长10.5mm、全长120mm
Angled, 10.5mm in blade length, 120mm in total length



巩膜隧道刀 Scleral Discission Knife



52277A

月形、刀宽2mm、全长113mm
Tunnel knife, 2mm in width of knife, 113mm in total length

裂隙穿刺刀 Slit Knife



52278A

角形、刀宽3mm、全长113mm
Angled, 3mm in width of knife, 113mm in total length

晶体植入切开刀 Discission Knife for IOL Implantation



52279A

角形、刀宽4mm、全长113mm
Angled, 4mm in blade width, 113mm in total length

角巩膜切开刀 Corneoscleral Knife



52286A

直形、斜角10°、全长111mm
Straight, 10°, 111mm in total length

巩膜穿刺刀 Scleral Discission Knife



52300

直形、刃长5mm、全长132mm
Straight, 5mm in blade length, 132mm in total length

劈核刀 Cystotome Knife



52283D


左角形、头宽0.7mm、全长110mm
Left angled, 0.7mm in tip width, 110mm in total length



52284D

右角形、头宽0.7mm、全长110mm
Right angled, 0.7mm in tip width, 110mm in total length

劈核刀 Cystotome Knife



52285D

圆头、刃长0.5mm、全长107mm
Round tip, 0.5mm in blade length, 107mm in total length

劈核刀 Cystotome Knife



52287D

左楔形刃长0.9mm、左式、全长110mm
Cuniform, 0.9mm in blade length, left, 110mm in total length

劈核刀 Cystotome Knife



52294D

刃长1.6mm、全长110mm
1.6mm in blade length, 110mm in total length



52295D

刃长1.8mm、全长110mm
1.8mm in blade length, 110mm in total length

虹膜复位器 Iris Spatula



52390

直形、头宽1.7mm、全长130mm
Straight, 1.7mm in tip width, 130mm in total length

52391

弯形、头宽1.7mm、全长130mm
Curved, 1.7mm in tip width, 130mm in total length

虹膜复位器 Iris Spatula



52429T

一头宽0.4mm、一头宽1.2mm
0.8mm in one tip width, 1.2mm in another tip width

虹膜复位分离器 Iris Spatula



52404A

短角形、头长11mm、全长120mm
Short angled, 11mm in tip length, 120mm in total length

52404T

短角形、头长11mm、全长120mm
Short angled, 11mm in tip length, 120mm in total length

虹膜复位器 Iris Spatula



52414

直形、头宽1.7mm、全长117mm
Straight, 1.7mm in tip width, 117mm in total length

52415

直形、头宽1.2mm、全长117mm
Straight, 1.2mm in tip width, 117mm in total length

虹膜复位分离器 Iris Spatula



52417A

长角形，头宽 $\phi 0.2\text{mm}$ 、全长112mm
Long angled, R45, $\phi 0.2\text{mm}$ in tip width, 112mm in total length

虹膜复位器（注水式） Iris Spatula (Infusion)



56135A

头部 $0.55 \times 0.8\text{mm}$
 $0.55 \times 0.8\text{mm}$ in tip length

56137A

头长7mm, $\phi 0.5\text{mm}$
7mm in tip length, $\phi 0.5\text{mm}$

白内障匙 Lens Spoon

Lens Spoon



52574A

匙长 $8 \times 2.5\text{mm}$ 、全长120mm
 $8 \times 2.5\text{mm}$ in spoon length, 120mm in total length

睑板腺囊肿锐匙（椭圆形） Chalazion Curette



52630

匙形 $2 \times 1\text{mm}$ 、全长130mm
Oval-shaped $2 \times 1\text{mm}$, 130mm in total length

52632

匙形 $3 \times 1.5\text{mm}$ 、全长130mm
Oval-shaped $3 \times 1.5\text{mm}$, 130mm in total length

52634

匙形 $3.5 \times 2\text{mm}$ 、全长130mm
Oval-shaped $3.5 \times 2\text{mm}$, 130mm in total length

52636

匙形 $4 \times 2.5\text{mm}$ 、全长130mm
Oval-shaped $4 \times 2.5\text{mm}$, 130mm in total length

52638

匙形 $5 \times 3\text{mm}$ 、全长130mm
Oval-shaped $5 \times 3\text{mm}$, 130mm in total length

52640

匙形 $6 \times 3.5\text{mm}$ 、全长130mm
Oval-shaped $6 \times 3.5\text{mm}$, 130mm in total length

睑板腺囊肿锐匙（圆形） Chalazion Curette (Round)



52650

匙径1mm、全长130mm
Curette $\phi 1\text{mm}$, 130mm in total length

52652

匙径1.5mm、全长130mm
Curette $\phi 1.5\text{mm}$, 130mm in total length

52653

匙径2mm、全长130mm
Curette $\phi 2\text{mm}$, 130mm in total length

52654

匙径2.5mm、全长130mm
Curette $\phi 2.5\text{mm}$, 130mm in total length

52655

匙径3mm、全长130mm
Curette $\phi 3\text{mm}$, 130mm in total length

52656

匙径3.5mm、全长130mm
Curette $\phi 3.5\text{mm}$, 130mm in total length

晶体线环 Lens Loop



52488A

弯形、头部4×6mm、全长120mm
Curved loop, tip 4mm (width)×6mm (length) curved loop,
120mm in total length



52489A

弯形有齿、头部3×7.5mm (倒齿)、全长122mm
Curved serrated, tip 3mm (width)×7.5mm (length),
122mm in total length



52490A

头部3×7.5mm
3×7.5mm tip

晶体线环 (注水式) Lens Loop (Infusion)



56133A

头部4×6mm
Tip 4mm (width)×6mm (length)

晶体线环 (注水式、有齿) Lens Loop (Infusion, serrated)



56136A

头部 (3.5×11)
3.5mm (width)×11mm (length) slightly

56167A

头部5×7mm
Tip 5mm (width)×7mm (length)

56168A

头部3×9mm
Tip 3mm (width)×9mm (length)

折叠式人工晶体植入器 Foldable IOL Implantation Holder



55171T

头部直径 φ 1.3mm
Tip φ 1.3mm

一体式人工晶体推注器 One-piece IOL Implantation Holder



55163T

头部1.7mm×0.7mm
Tip 1.7mm×0.7mm

人工晶体植入镊 IOL Implantation Forceps



53316A

角形、头部直径 φ 0.6mm、50°、全长105mm
Tip φ 0.6mm, angled 50°, 105mm in total length

人工晶体植入镊 IOL Implantation Forceps



53316T

角形、头部直径 φ 0.6mm、50°、全长114mm
Tip φ 0.6mm, angled 50°, 114mm in total length

人工晶体植入镊（折叠式） IOL Implantation Forceps (Foldable)



53332T

夹持面长6.5mm、30°、全长108mm，带锁止扣
6.5mm in length of long platform, angled 30°,
108mm in total length

人工晶体植入镊 IOL Implantation Forceps



53580A

头部弯形、R12、全长112mm

53580T

Curved, R12, 112mm in total length
头部弯形、R12、全长112mm
Curved, R12, 112mm in total length

人工晶体植入镊 IOL Implantation Forceps



53591T

圆板形、全长112mm
Round tip, 112mm in total length



53592T

L形、全长112mm
L-shaped, 112mm in total length

撕囊镊 Capsulorhexis Forceps



53317A

弧形圆头、宽0.3mm、全长105mm
Curved, round tips, 0.3mm in tip width,
105mm in total length



53383A

弧形角形头、头部宽0.5mm、全长105mm
Curved, angled, 0.5mm in tip width,
105mm in total length

撕囊镊 Capsulorhexis Forceps



53384T

直形圆角头、锐钩宽0.65mm、全长105mm
Straight, in the width of sharp hook,
105mm in total length

撕囊镊 Capsulorhexis Forceps



53589T

弧形角形头、头部宽0.5mm、全长114mm
Curved, angled, 0.5mm in tip width, 114mm in total length

撕囊镊 Capsulorhexis Forceps



53381

弧形瓜子头、头宽0.36mm、防滑线、全长110mm
Curved, claw-shaped tip, 0.36mm in tip width,
110mm in total length

撕囊镊 Capsulorhexis Forceps



53382

弧形角形头（弹簧式）、头宽0.6mm、防滑线、全长118mm
Curved, angled, spring type, 0.6mm in tip width, 118mm in total length

微创撕囊镊 Small Incision Copsulorhexis Forceps



53177TS

弯角形头，全长155mm
Curved, angled tip, 155mm in total length



旋转式微创撕囊镊 Rotation Small Incision Copsulorhexis Forceps



53179TX

头宽0.15mm，全长142mm
0.15mm in head width, 142mm in total length



缝线结扎镊 Suturing Forceps



53325A

直平台长4.5mm、唇头钩0.15mm、全长105mm
Straight, platform 4.5mm in length, 1×2 teeth of 0.15mm, 105mm in total length

53326A

直平台长4.5mm、唇头钩0.2mm、全长105mm
Straight, platform 4.5mm in length, 1×2 teeth of 0.2mm, 105mm in total length

53329A

直平台长4.5mm、唇头钩0.12mm、全长105mm
Straight, platform 4.5mm in length, 1×2 teeth of 0.12mm, 105mm in total length



缝线结扎镊 Suturing Forceps



53325T

直平台长4.5mm、唇头钩0.15mm、全长114mm
Straight, platform 4.5mm in length, 1×2 teeth of 0.15mm, 114mm in total length

53326T

直平台长4.5mm、唇头钩0.2mm、全长114mm
Straight, platform 4.5mm in length, 1×2 teeth of 0.2mm, 114mm in total length

53329T

直平台长4.5mm、唇头钩0.12mm、全长114mm
Straight, platform 4.5mm in length, 1×2 teeth of 0.12mm, 114mm in total length



缝线结扎镊 Suturing Forceps



53411A

直平台长4.5mm、0.15mm唇头钩、全长114mm
Straight, platform 4.5mm in length, 1×2 teeth of 0.15mm, 114mm in total length



53412A

直平台长4.5mm、0.20mm唇头钩、全长114mm
Straight, platform 4.5mm in length, 1×2 teeth of 0.20mm, 114mm in total length



53413A

直平台长4.5mm、0.12mm唇头钩、全长114mm
Straight, platform 4.5mm in length, 1×2 teeth of 0.12mm, 114mm in total length



缝线结扎镊 Suturing Forceps



53570A

直平台长6.5mm、0.4mm唇头钩、全长105mm
Straight, platform 6.5mm in length, 1×2 teeth of 0.4mm, 105mm in total length



缝线结扎镊 Suturing Forceps



53574A

弹簧式、直平台长5.0mm、0.15mm唇头钩、全长110mm
Spring type, straight, platform 5.0mm in length, 1×2 teeth of 0.15mm, 110mm in total length



缝线结扎镊 Suturing Forceps



53583A

53584A



直平台长4.5mm、0.15mm唇头钩、全长88mm
Straight, platform 4.5mm in length, 1×2 teeth of 0.15mm, 88mm in total length

缝线结扎镊 Suturing Forceps



53593T

53583T

53584T



直平台长4.5mm、0.12mm唇头钩、全长88mm
Straight, platform 4.5mm in length, 1×2 teeth of 0.12mm, 88mm in total length

直平台长4.5mm、0.15mm唇头钩、全长88mm
Straight, platform 4.5mm in length, 1×2 teeth of 0.15mm, 88mm in total length

直平台长4.5mm、0.20mm唇头钩、全长88mm
Straight, platform 4.5mm in length, 1×2 teeth of 0.20mm, 88mm in total length

系结镊 Tying Forceps



53320A



直平台长4.5mm、全长105mm
Straight, platform 4.5mm in length, 105mm in total length

系结镊 Tying Forceps



53320T



直平台长4.5mm、全长114mm
Straight, platform 4.5mm in length, 114mm in total length

系结镊 Tying Forceps



53322A



弯平台长4.0mm、全长105mm
Curved, platform 4.0mm in length, 105mm in total length

系结镊 Tying Forceps



53322T



弯平台长4.0mm、全长114mm
Curved, platform 4.0mm in length, 114mm in total length

系结镊 Tying Forceps



53323A

53324A



长角平台长4.5mm、全长105mm
Angled, platform 4.5mm in length, 105mm in total length

短角平台长4.5mm、全长114mm
Angled, platform 4.5mm in length, 114mm in total length

系结镊 Tying Forceps



53323T

53324T



长角平台长4.5mm、全长105mm
Angled, platform 4.5mm in length, 105mm in total length

短角平台长4.5mm、全长114mm
Angled, platform 4.5mm in length, 114mm in total length

系结镊 Tying Forceps



53414A

直平台长4.5mm、全长114mm
Straight, platform 4.5mm in length, 114mm in total length



系结镊 Tying Forceps



53415A

弯平台长4.5mm、全长114mm
Curved, platform 4.5mm in length, 114mm in total length



系结镊 Tying Forceps



53416A

长角平台长4.5mm、全长114mm
Angled platform 4.5mm in length, 114mm in total length



53417A

短角平台长4.5mm、全长114mm
Angled, platform 4.5mm in length, 114mm in total length



系结镊 Tying Forceps



53576A

弹簧式、直平台长5.0mm、全长110mm
Spring type, straight, platform 5.0mm in length, 110mm in total length



系结镊 Corneal Tying Forceps



53586A

直平台长4.5mm、全长88mm
Straight, platform 4.5mm in length, 88mm in total length



53586T

直平台长4.5mm、全长88mm
Straight, platform 4.5mm in length, 88mm in total length



无损伤镊 Corneal Scleral Forceps



53288A

直半圆、全长105mm
Straight, semicircle tips, 105mm in total length



无损伤镊 Corneal Scleral Forceps



53288T

直半圆、全长114mm
Straight shafts with micro pierce type tips, 114mm in total length



无损伤镊 Corneal Forceps



53585A

直半圆、全长88mm
Straight, semicircle tips, 88mm in total length



53585T

直半圆、全长88mm
Straight, semicircle tips, 88mm in total length



虹膜镊 Iris Forceps



53004

曲1×2唇头钩、全长70mm
Bend 1×2 teeth, 70mm in total length



53042

弯1×2唇头钩、全长60mm
Curved 1×2 teeth, 60mm in total length

眼用镊 Ophthalmic Forceps



53070

直唇头齿、齿距0.5mm、全长100mm
Straight, 0.5mm in serrated tips distance, 100mm in total length



53072

弯唇头齿、齿距0.5mm、全长100mm
Curved, 0.5mm in serrated tips distance, 100mm in total length



53076

直、1×2唇头钩、全长100mm
Straight, 1×2 teeth, 100mm in total length



53078

弯、1×2唇头钩、全长100mm
Curved, 1×2 teeth, 100mm in total length

结膜镊 Conjunctival Forceps



53120

直唇头齿、齿距0.4mm、全长70mm
Straight, 0.4mm in serrated tips distance, 70mm in total length



53124

直1×2唇头钩、全长70mm
Straight, 1×2 teeth, 70mm in total length

结膜镊 Conjunctival Forceps



53132

弯唇头齿、全长70mm
Curved, 70mm in total length



53138

弯1×2唇头齿、全长70mm
Curved, 1×2 teeth, 70mm in total length

异物镊 Foreign Body Forceps



53170

无唇头齿、全长100mm
Sharp pointed tips, 100mm in total length

晶体囊镊 Capsulotomy Forceps



53274

弯长方匙、匙径1.2×3mm、全长70mm
Oblong spoon tips, curette 1.2×3mm, 70mm in total length

固定镊 Fixation Forceps



53340

直形、齿距0.7mm、全长110mm
Straight, 0.7mm in serrated tips distance, 110mm in total length

睑板腺囊肿镊 Chalazion Forceps



53430

椭圆形、20×13mm、全长90mm
Oval blades, 20×13mm, 90mm in total length

53432

椭圆形、25×16mm、全长90mm
Oval blades, 25×16mm, 90mm in total length

睑板腺囊肿镊 Chalazion Forceps



53440

圆环形、Φ12mm、全长90mm
Round blades, Φ12mm, 90mm in total length

睫毛镊 Cilia Forceps



53530

圆头、夹持面长1.5mm、全长90mm
Round tips, long platform 1.5mm in length, 90mm in total length

显微眼用镊 Micro-ophthalmic Forceps



53327A

直唇头齿、齿距0.25mm、全长105mm
Straight, serrated tips, 0.25mm in distance, 105mm in total length

53327T

直唇头齿、齿距0.25mm、全长114mm
Straight, serrated tips, 0.25mm in distance, 114mm in total length

移核镊 Capsul Extractive Forceps



53328A

角形、唇头齿长8mm、全长105mm
Angled, 8mm in serrated length, 105mm in total length

巩膜塞夹持镊 Sclera Plug Holding Forceps



53582A

反力式、夹持槽长2mm、全长110mm
Curved, grasping trough 2mm in length, 110mm in total length

结膜剪 Conjunctival scissors



54095A

弯形 刃长17mm 全长115mm
Curved, 17mm in blade length, 115mm in total length

角膜剪 Corneal Scissors



54102

弯形 R5.5mm 全长90mm
Curved, R5.5mm, 90mm in total length

角膜剪 Corneal Scissors



54053B

直形、刃长14mm、全长120mm
Straight, 14mm in blade length, 120mm in total length

小梁剪 Micro-corneal Scissors



54108A

弯尖、刃长12mm、全长120mm
Curved, 12mm in blade length,
120mm in total length

54108B

弯尖、刃长12mm、全长120mm
Curved, 12mm in blade length,
120mm in total length

54110A

直尖、刃长12mm、全长120mm
Straight tips, 12mm in blade length,
120mm in total length

角膜剪 Corneal Scissors



54109A

弯形、刃长14mm、全长120mm
Slight curved, 14mm in blade length,
120mm in total length

54109B

弯形、刃长14mm、全长120mm
Slight curved, 14mm in blade length,
120mm in total length

膜状内障剪 Cataract Scissors Membrane



54264A

直形、刃长13mm、全长117mm
Straight tips, 13mm in blade length, 117mm in total length

54264B

直形、刃长13mm、全长117mm
Straight tips, 13mm in blade length, 117mm in total length

膜状内障剪 Cataract Scissors Membrane



54268A

直形、刃长16mm、全长120mm
Straight, 16mm in blade length, 120mm in total length

54268B

直形、刃长16mm、全长120mm
Straight, 16mm in blade length, 120mm in total length

膜状内障剪 Cataract Scissors Membrane



54269A

弯形、刃长16mm、全长120mm
Curved, 16mm in blade length, 120mm in total length

54269B

弯形、刃长16mm、全长120mm
Curved, 16mm in blade length, 120mm in total length

维纳斯剪 Vannas Scissors



54136B

直形、刃长6.5mm、全长84mm
Straight, 6.5mm in blade length, 84mm in total length

54137B

尖头、刃长8mm、全长84mm
Sharp tips, 8mm in blade length, 84mm in total length

54138B

直形、刃长6.5mm、全长84mm
Straight, 6.5mm in blade length, 84mm in total length

维纳斯剪 Vannas Scissors



54139B

弯形、刃长8mm、全长90mm
Curved, 8mm in blade length, 90mm in total length



54140B

弯形、刃长6.5mm、全长84mm
Curved, 6.5mm in blade length, 84mm in total length



54156B

弯形、刃长5mm、全长80mm
Curved, 5mm in blade length, 80mm in total length



54050B

直形、刃长8mm、全长84mm
Straight, 8mm in blade length, 84mm in total length

54147B

弯形、刃长9mm、全长84mm
Curved, 9mm in blade length, 84mm in total length



眼用手术剪 Ophthalmic Scissors



54002

直尖头、全长100mm
Straight sharp pointed tips, 100mm in total length

54006

直圆头、全长100mm
Straight blunt tips, 100mm in total length

54010

弯尖头、全长100mm
Curved sharp pointed tips, 100mm in total length

54014

弯圆头、全长100mm
Curved blunt tips, 100mm in total length

眼用显微手术剪 Micro Ophthalmic Scissors



54024

弯尖头、全长85mm
Curved sharp pointed tips, 85mm in total length

54026

弯圆头、全长85mm
Curved round pointed tips, 85mm in total length

54020

直尖头、全长85mm
Straight sharp pointed tips, 85mm in total length

54022

直圆头、全长85mm
Curved round pointed tips, 85mm in total length

虹膜剪 Iris Scissors



54060

双尖、全长90mm
Double pointed tips, 90mm in total length

显微持针钳 Needle Scissors



55115

直、小方柄、全长130mm
Straight, 130mm in total length

显微持针钳 Micro-needle Holder



55152B

弯形、全长110mm
Curved, 110mm in total length

55152T

弯形、全长110mm
Curved, 110mm in total length

显微持针钳 Micro-needle Holder



55153T

弯形、全长118mm
Curved, 118mm in total length

显微持针钳 Micro-needle Holder



55154A



直狭头、全长110mm
Straight, 110mm in total length

显微持针钳 Micro-needle Holder



55158T



弯形、全长110mm
Curved, 110mm in total length

显微持针钳 Micro-needle Holder



55159A



弯锁止钩、全长140mm
Curved, with lock, 140mm in total length

显微持针钳 Micro-needle Holder



55161T



弯锁止钩、全长110mm
Round curved handle, with lock, 110mm in total length

白内障注吸器 Irrigation/aspiration Cannula

56125A



弯双管 $\phi 0.6\text{mm}$
 $\phi 0.6\text{mm}$, with adapter, curved double tubes

12点位置注吸器 12 O' clock Cannula



56183

左式 $\phi 0.5\text{mm}$
Left $\phi 0.5\text{mm}$

56184

右式 $\phi 0.5\text{mm}$
Right $\phi 0.5\text{mm}$



白内障注吸器 Coaxial I/A Cannula



56188A

弯三接式, $\phi 0.6\text{mm}$
Curved, triple joint, $\phi 0.6\text{mm}$



白内障注吸器 Coaxial I/A Cannula



56190

直、三接式、 $\phi 0.6\text{mm}$
Straight tip, triple joint, $\phi 0.6\text{mm}$



白内障注吸器 Coaxial I/A Cannula



56191A(含56199A)

(Conclude 56199A)

弯双管 $\phi 0.7\text{mm}$
Double curved tubes, $\phi 0.7\text{mm}$



12点位置注吸器

12 O' clock Cannula



56195

左式 $\phi 0.7\text{mm}$
Left $\phi 0.7\text{mm}$

56196

右式 $\phi 0.7\text{mm}$
Right $\phi 0.7\text{mm}$

白内障注吸器

Irrigation/aspiration Units



56216A

弯双管 $\phi 0.6\text{mm}$

Double curved tubes $\phi 0.6\text{mm}$

56217A

弯双管 $\phi 0.8\text{mm}$

Double curved tubes $\phi 0.8\text{mm}$

56218A

弯双管 $\phi 0.9\text{mm}$

Double curved tubes $\phi 0.9\text{mm}$

冲洗针（弯泪道注水式）

Aspiration Needle



56140

$\phi 0.45\text{mm}$ $\phi 0.6\text{mm}$ $\phi 0.8\text{mm}$

冲洗针（直泪道注水式）

Irrigation Needle



56141

$\phi 0.6\text{mm}$ $\phi 0.7\text{mm}$ $\phi 0.8\text{mm}$
 $\phi 0.9\text{mm}$ $\phi 1.0\text{mm}$ $\phi 1.1\text{mm}$

冲洗针（弯泪道注水式）

Aspiration Needle Curved



56142

$\phi 0.6\text{mm}$ $\phi 0.7\text{mm}$ $\phi 0.8\text{mm}$ $\phi 0.9\text{mm}$
 $\phi 1.0\text{mm}$ $\phi 1.1\text{mm}$ $\phi 1.2\text{mm}$

冲洗针

Aspiration Needle



56144

直形 $\phi 0.4\text{mm}$ 全长27mm
Straight $\phi 0.4\text{mm}$ tip 27mm in total length



56146

角形 $\phi 0.4\text{mm}$ 15° 全长27mm
Angled $\phi 0.4\text{mm}$ tip 27mm in total length



56151A

角形30° $\phi 0.3\text{mm}$ 全长18mm
Angled 30° $\phi 0.3\text{mm}$ tip 18mm in total length



56152A

角形60° $\phi 0.3\text{mm}$ 全长15mm
Angled 60° $\phi 0.3\text{mm}$ tip 15mm in total length

冲洗针

Aspiration Needle



56149A

弯形 $\phi 0.7\text{mm}$ 全长26mm
Curved $\phi 0.7\text{mm}$, 26mm in total length

注水器 Aspirating Needle



56197A



角形 头宽1.1×0.5mm
Angled, 1.1×0.5mm tip

角膜上皮剥离器 Corneal Spatula



52423A

圆头, 头宽2.5mm、全长118mm
Round, 2.5mm in head width, 118mm in total length



52424A

半圆头宽1.5mm、全长112mm
Semicircular head, 1.5mm in head width, 112mm in total length

角膜上皮钩 Corneal Epidermic Hook



52437A



头宽1mm、全长120mm
1mm in head width, 120mm in total length

角膜上皮铲 Corneal Epidermic Pick



52621A

T形、头宽1.8mm、全长116mm
T-shaped, 1.8mm in head width, 116mm in total length

52624A

方头、头宽1.8mm、全长116mm
Square head, 1.8mm in head width, 116mm in total length

角膜上皮镊 Corneal Epidermic Forceps



53387A



圆形 $\phi 0.6\text{mm}$ 、全长114mm
Round, $\phi 0.6\text{mm}$, 114mm in total length

角膜上皮镊 Corneal Epidermic Forceps



53595T



环形头 $\phi 1.1\text{mm}$ 、全长114mm
 $\phi 1.1\text{mm}$ circle, 114mm in total length

角膜上皮环钻 Corneal Epidermic Loop



52482B



$\phi 7.5\text{mm}$ 、全长115mm
 $\phi 7.5\text{mm}$, 115mm in total length

52483B

$\phi 8.0\text{mm}$ 、全长115mm
 $\phi 8.0\text{mm}$, 115mm in total length

52484B

$\phi 8.5\text{mm}$ 、全长115mm
 $\phi 8.5\text{mm}$, 115mm in total length

52485B

$\phi 9.0\text{mm}$ 、全长115mm
 $\phi 9.0\text{mm}$, 115mm in total length

角膜上皮酒精环 Corneal Epidermic Loop



52502B



内径 $\phi 8.0\text{mm}$ 、全长115mm
 $\phi 8.0\text{mm}$, 115mm in total length

52503B

内径 $\phi 8.5\text{mm}$ 、全长115mm
 $\phi 8.5\text{mm}$, 115mm in total length

52504B

内径 $\phi 9.0\text{mm}$ 、全长115mm
 $\phi 9.0\text{mm}$, 115mm in total length

52505B

内径 $\phi 9.5\text{mm}$ 、全长115mm
 $\phi 9.5\text{mm}$, 115mm in total length

角膜环钻(有芯) Trephines



56003

刃部口径3mm、全长62mm
Diameter of blade part 3mm, 62mm in total length

56004

刃部口径3.5mm、全长62mm
Diameter of blade part 3.5mm, 62mm in total length

56005

刃部口径4mm、全长62mm
Diameter of blade part 4mm, 62mm in total length

56006

刃部口径4.5mm、全长62mm
Diameter of blade part 4.5mm, 62mm in total length

56007

刃部口径5mm、全长62mm
Diameter of blade part 5mm, 62mm in total length

56008

刃部口径5.5mm、全长62mm
Diameter of blade part 5.5mm, 62mm in total length

56009

刃部口径6mm、全长62mm
Diameter of blade part 6mm, 62mm in total length

56010

刃部口径6.5mm、全长62mm
Diameter of blade part 6.5mm, 62mm in total length

56011

刃部口径7mm、全长62mm
Diameter of blade part 7mm, 62mm in total length

56012

刃部口径8mm、全长62mm
Diameter of blade part 8mm, 62mm in total length



角膜环钻(无芯) Trephines



56050A

刃部口径6.0mm、全长60mm
Diameter of blade part 6.0mm, 60mm in total length

56051A

刃部口径6.25mm、全长60mm
Diameter of blade part 6.25mm, 60mm in total length

56052A

刃部口径6.5mm、全长60mm
Diameter of blade part 6.5mm, 60mm in total length

56053A

刃部口径6.75mm、全长60mm
Diameter of blade part 6.75mm, 60mm in total length

56054A

刃部口径7mm、全长60mm
Diameter of blade part 7mm, 60mm in total length

56055A

刃部口径7.25mm、全长60mm
Diameter of blade part 7.25mm, 60mm in total length

56056A

刃部口径7.5mm、全长60mm
Diameter of blade part 7.5mm, 60mm in total length

56057A

刃部口径7.75mm、全长60mm
Diameter of blade part 7.75mm, 60mm in total length

56058A

刃部口径8mm、全长60mm
Diameter of blade part 8mm, 60mm in total length

56059A

刃部口径8.25mm、全长60mm
Diameter of blade part 8.25mm, 60mm in total length



角膜环钻(无芯) Trephines



56306

刃部口径3mm、全长22mm
Diameter of blade part 3mm, 22mm in total length

56307

刃部口径3.25mm、全长22mm
Diameter of blade part 3.25mm, 22mm in total length

56308

刃部口径3.5mm、全长22mm
Diameter of blade part 3.5mm, 22mm in total length

56309

刃部口径3.75mm、全长22mm
Diameter of blade part 3.75mm, 22mm in total length

56310

刃部口径4mm、全长22mm
Diameter of blade part 4mm, 22mm in total length

56311

刃部口径4.25mm、全长26mm
Diameter of blade part 4.25mm, 26mm in total length

56312

刃部口径4.5mm、全长26mm
Diameter of blade part 4.5mm, 26mm in total length

56313

刃部口径4.75mm、全长26mm
Diameter of blade part 4.75mm, 26mm in total length

56314

刃部口径5mm、全长26mm
Diameter of blade part 5mm, 26mm in total length

56315

刃部口径5.25mm、全长26mm
Diameter of blade part 5.25mm, 26mm in total length

56316

刃部口径5.5mm、全长26mm
Diameter of blade part 5.5mm, 26mm in total length

56317

刃部口径5.75mm、全长26mm
Diameter of blade part 5.75mm, 26mm in total length

56318

刃部口径6mm、全长26mm
Diameter of blade part 6mm, 26mm in total length

56319

刃部口径6.25mm、全长26mm
Diameter of blade part 6.25mm, 26mm in total length

56320

刃部口径6.5mm、全长26mm
Diameter of blade part 6.5mm, 26mm in total length

56321

刃部口径6.75mm、全长26mm
Diameter of blade part 6.75mm, 26mm in total length

56322

刃部口径7mm、全长26mm
Diameter of blade part 7mm, 26mm in total length

56323

刃部口径7.25mm、全长30mm
Diameter of blade part 7.25mm, 30mm in total length

56324

刃部口径7.5mm、全长30mm
Diameter of blade part 7.5mm, 30mm in total length

56325

刃部口径7.75mm、全长30mm
Diameter of blade part 7.75mm, 30mm in total length

56326

刃部口径8mm、全长30mm
Diameter of blade part 8mm, 30mm in total length

56327

刃部口径8.25mm、全长30mm
Diameter of blade part 8.25mm, 30mm in total length

56328

刃部口径8.5mm、全长30mm
Diameter of blade part 8.5mm, 30mm in total length

泪囊牵开器 Lachrymal Sac Retractor



51430

3×3爪、全长82mm
3×3 claws, 82mm in total length

核抛光器 Capsule Polisher



51472A

角形头部0.6×0.8mm、全长113mm
Angled, 0.6×0.8 mm elliptical tip, 113mm in total length

晶体切核器 Lens Nucleus Manipulator



51473A

角形头宽0.25mm、全长112mm
Angled, 0.25mm in tip width, 112mm in total length

异物针 Foreign Body Needle



52260

直形、刃长5mm、全长125mm
Straight, 5mm in blade length, 125mm in total length



52262

弯形、刃长5mm、全长125mm
Curved, 5mm in blade length, 125mm in total length

碎核托板 Nucleus Spatula



52334

左式、头部宽2.8mm、全长125mm
Left, 2.8mm in head width, 125mm in total length



52335

右式、头部宽2.8mm、全长125mm
Right, 2.8mm in head width, 125mm in total length

泪点扩张器 Lachrymal Dilator



56020

短锥、全长110mm
Short taper, 110mm in total length



56022

长锥、全长110mm
Long taper, 110mm in total length

泪道探针 Lachrymal Probe



56045

φ0.5-0.6mm 全长130mm
130mm in total length

材质:银合金
Material: Silver alloy

56046

φ0.7-0.8mm 全长130mm
130mm in total length

材质:银合金
Material: Silver alloy

56047

φ0.9-1.0mm 全长130mm
130mm in total length

材质:银合金
Material: Silver alloy

56048

φ1.1-1.2mm 全长130mm
130mm in total length

材质:银合金
Material: Silver alloy

小梁切开刀 Trabeculotomy Probe



56090

左式 φ0.2mm、全长57mm
Left φ0.2mm, 57mm in total length

56092

右式 φ0.2mm、全长57mm
Right φ0.2mm, 57mm in total length

巩膜咬切器 Corneoscleral Punch



56100

刃宽1.5mm 全长85mm
1.5mm in blade width, 85mm in total length

56102

刃宽2mm 全长85mm
2mm in blade width, 85mm in total length

56104

刃宽3mm 全长85mm
3mm in blade width, 85mm in total length

小梁咬切器 Micro trabecula punche



56105A
56106A

φ 1mm、镊式半圆
φ 1mm,semicircle,forceps-type
φ 1.5mm、镊式半圆
φ 1.5mm,semicircle,forceps-type



蝶形注液针 Infusion Cannulas



56156
56157
56158
56159
56160

头长2.5mm,全长8.1mm
2.5mm in beveled tips length,
8.1mm in total length
头长3mm,全长8.6mm
3mm in beveled tips length,
8.6mm in total length
头长4mm,全长9.6mm
4mm in beveled tips length,
9.6mm in total length
头长5mm,全长10mm
5mm in beveled tips length,
10mm in total length
头长6mm,全长11.6mm
6mm in beveled tips length,
11.6mm in total length



烧灼止血器 Cautery

52700

球针型、圆形 φ 5.5mm、全长119mm
Ball shape tip, φ 5.5mm, 119mm in total length



笛针 Aspiraion Needle



56170A

直形、头长35mm、全长140mm
Straight, 35mm in tip length, 140mm in total length



刀片夹持器 Blade Holder



56202B
56202T

弹簧式、全长120mm
With spring, 120mm in total length

弹簧式、全长120mm
With spring, 120mm in total length



眼睑板 Eyelid Plate



51494

全长100mm
100mm in total length

巩膜塞 Scleral Plug

56500

φ 0.9mm
φ 0.9mm

全长3.5mm
3.5mm in total length



注水手柄 Infusion Handpiece

56203

双接头 全长100mm
100mm in total length

56205

单锥头 全长50mm
50mm in total length



注水手柄 Infusion Handpiece of Cataract

56206T

单锥头 全长100mm
Single head of cone, 100mm in total length



眼用规 Castroviejo Caliper



56210

固定式 规格3.5/4.5/5.0/6.0mm
Castroviejo caliper 3.5/4.5/5.0/6.0mm

眼用规 Castroviejo Caliper



56211A

可调试测量范围0~20
Adjustable, measuring from 0~20 at 0.1mm increment

眼用规 Castroviejo Caliper



56215

固定式、规格3.5/4mm
Fixation, 3.5mm/4mm tips

56215A

固定式、规格5/5.54mm
Fixation, 5mm/5.5mm tips

消毒盒 Sterilizing Box




57020A

不锈钢、205×125×34mm
Stainless steel, 205×125×34mm

57021A

不锈钢、235×150×34mm
Stainless steel, 235×150×34mm

SYX6A显微手术镊器械包 SYX6A Micro Surgical Forceps Set




序号	产品名称	代号	形式规格
1	系结镊 Tying Forceps	53323A	角平台长4.5mm,全长105mm Long angled, platform 4.5mm in length, 105mm in total length
2	无损伤镊 Corneal Scleral Forceps	53288A	直半圆,全长105mm Straight shafts with micro pierce type tips, 105mm in total length
3	缝线结扎镊 Suturing Forceps	53325A	直平台长4.5mm、0.15mm唇头钩,全长105mm Straight, platform 4.5mm in length, 1×2 teeth of 0.15mm 105mm in total length
4	显微眼用镊 Micro-ophthalmic Forceps	53327A	直唇头齿、齿距0.25mm,全长105mm Straight, serrated tips, 0.25mm in distance, 105mm in total length
5	系结镊 Tying Forceps	53320A	直平台长4.5mm、全长105mm Straight, platform 4.5mm in length, 105mm in total length
6	缝线结扎镊 Suturing Forceps	53326A	直平台,长4.5mm,0.20唇头钩,全长105mm Straight, platform 4.5mm in length, 1×2 teeth of 0.20mm, 105mm in total length

序号	产品名称	代号	形式规格
1	显微持针钳 Micro-needle Holder	55152B	弯形,全长110mm Curved, 110mm in total length
2	小梁剪 Micro-corneal Scissors	54108A/B	弯形,刃长12mm,全长120mm Curved, 12mm in blade length, 120mm in total length
3	刀片夹持器 Blade Holder	56202B	弹簧式,全长120mm With spring, 120mm in total length
4	显微眼用镊 Micro-ophthalmic Forceps	53327A	直唇头齿齿距0.25mm,全长105mm Straight, serrated tips, 0.25mm in distance, 105mm in total length
5	虹膜复位器 Iris Spatula	52391	弯形头宽1.7mm,全长130mm Curved, 1.7mm in tip width, 130mm in total length
6	缝线结扎镊 Suturing Forceps	53325A	直平台长4.5mm,0.15mm唇头钩,全长105mm Straight, platform 4.5mm in length, 1×2 teeth of 0.15mm, 105mm in total length
7-8	系结镊 Tying Forceps	53322A	弯平台长4.0mm,全长105mm Curved, platform 4.0mm in length, 105mm in total length
		53320A	直平台长4.5mm,全长105mm Straight, platform 4.5mm in length, 105mm in total length

SYX8实用显微手术器械包 SYX8 Practical Micro-ophthalmic Surgical Instrument Set



SX7 飞秒激光手术器械包 Femtosecond laser surgery instruments set



序号	产品名称	代号	形式规格
1	开睑器 Eye Speculum	51288A	钢丝式,扩张杆长15mm,全长50mm Wire speculum, stainless 15mm blade in length, 50mm in total length
2	冲洗针 15° Aspiration Needle Angled 15°	56146	角形 $\phi 0.4\text{mm}$ 全长27mm $\phi 0.4\text{mm}$ tip, 27mm in total length
3	系结镊 Tying Forceps	53338DH	直形, 全长114mm Straight, 114mm in total length
4	虹膜复位器 Iris Spatula	52429T	双头式, 全长140mm, 一头0.8mm, 另一头1.2mm Double-headed, Head width of 0.8mm and 1.2mm, 140mm in total length
5	飞秒透镜分离器 Femtosecond laser surgery lens separator Double-headed	52395T	双头, 左端 $\phi 0.3\text{mm}$, 右端 $\phi 0.2\text{mm}$ $\phi 0.3\text{mm}$ and $\phi 0.2\text{mm}$
6	飞秒分离匙 Femtosecond laser surgery corneal flap separator	52435T	头端 $\phi 0.3\text{mm}$ $\phi 0.3\text{mm}$ tip
7	飞秒透镜镊 Femtosecond laser surgery lens Forcep	53155T	弯形 (网纹) Curved Textured

SX9 人工晶体植入显微手术器械包 SYX9 Micro-Operation Instruments Set for Intraocular Lens Implantation



序号	产品名称	代号	形式规格
1	膜状内障剪 Cataract Scissors Membrane	54264B	直形, 刃长10mm, 全长117mm Straight tips, 10mm in blade length, 117mm in total length
2	角膜剪 Corneal Scissors	54109B	弯形, 刃长14mm, 全长120mm Slight curved, 14mm in blade length, 120mm in total length
3	系结镊 Tying Forceps	53576A	弹簧式, 直平台长5.0mm, 全长110mm Spring type, straight, platform 5.0mm in length, 110mm in total length
4	缝线结扎镊 Suturing Forceps	53574A	弹簧式, 直平台长5.0mm, 0.15mm头钩, 全长110mm Straight, platform 5.0mm in length, 110mm in total length
5	系结镊 Tying Forceps	53323A	角平台, 平台长4.5mm, 全长105mm Long angled, platform 4.5mm in length, 105mm in total length
6	白内障注吸器 Coaxial I/A Cannula	56191A	弯双管, $\phi 0.7\text{mm}$ Double curved tubes, $\phi 0.7\text{mm}$
7	开睑器 Eye Speculum	51288A	钢丝式, 扩张杆长15mm, 全长50mm Wire speculum, stainless 15mm blade in length, 50mm in total length
8	人工晶体定位钩 IOL Fixation Hook	51478A	弯圆 $\phi 0.2\text{mm}$, 15° 全长110mm 0.2mm cylindrical tip, angled 15°, 110mm in total length
9	显微持针钳 Micro-needle Holder	55152A/B	弯形, 全长110mm Curved, 110mm in total length

SYX14白内障小切口手术器械包 SYX14 Cataract Small Cut Surgery Instruments Set



序号	产品名称	代号	形式规格
1	月形隧道刀 Moon-shaped Dicsision Knife	50102T	刀头宽2.8mm、厚0.2mm 2.8mm in tip width, 0.2mm in thickness
2	弯头穿刺刀 Puncture Knife	50403T	刀头宽3.0mm、厚0.2mm 3.0mm in tip width, 0.2mm in thickness
3	膜状内障剪 Cataract Scissor Membrane	54269B	弯形、刃长16mm、全长120mm Curved, 16mm in blade length, 120mm in total length
4	小梁剪 Micro Corneal Scissor	54108B	弯尖头、刃长12mm、全长120mm Curved round tip, 12mm in blade length, 120mm in total length
5	撕囊镊 Capsulorhexis Forcep	53383A	弧形角形头、头宽0.5mm、 全长105mm Curved and angled tip, 0.5mm in tip width, 105mm in total length
6	缝线结扎镊 Suturing Forcep	53325A	唇头钩0.15mm、全长105mm 1×2 teeth of 0.15mm, 105mm in total length
7	系结镊 Tying Forcep	53323A	长角形平台长4.5mm、全长105mm Long angled platform, 4.5mm in length, 105mm in total length
8	碎核垫圈 Penetration Ring	52490A	头部3×7.5mm 3×7.5mm in head
9	碎核刀 Penetration Knife	52276D	角形、刃长10.5mm、全长120mm Angled, 10.5mm in blade length, 120mm in total length
10	注水手柄 Infusion Handle	56203	双接头、全长100mm Double connector, 100mm in total length

序号	产品名称	代号	形式规格
11	白内障注吸器 Cataract Irrigation/ aspiration Cannula	56191A	弯双管、Φ0.7mm Curved double cannula, Φ0.7mm
12	晶体线环 Lens Loop	56136A	注水式、头部3.5×11mm、 全长26mm Infusion, 3.5×11mm tip, 26mm in total length
13	开睑器 Speculum	51288A	钢丝式、扩张杆长15mm、 全长50mm Wire Speculum, stainless 15mm blade in length, 50mm in total length
14	眼用规 Eye Gauge	56210	固定式、圆形II Fixation, round II

SYX15青光眼显微手术器械包 SYX15 Glaucoma Micro-operation Instruments Set




序号	产品名称	代号	形式规格
1	冲洗针 Aspiration Needle	56146	弯形、Φ0.4mm、全长27mm Angled Φ0.4mm, tip 27mm in total length
2	冲洗针 Aspiration Needle	56144	直形、Φ0.4mm、全长27mm Straight Φ0.4mm, tip 27mm in total length
3	前房角切开刀 Goniotomy Knife	52231A	直形、刃长2mm、全长115mm Straight, 2mm in blade length, 115mm in total length
4	虹膜复位器 Iris Spatula	52415	直形、头宽1.2mm、全长117mm Straight, 1.2mm in tip width, 117mm in total length
5	虹膜复位器 Iris Spatula	52414	直形、头宽1.7mm、全长117mm Straight, 1.7mm in tip width, 117mm in total length
6	虹膜复位分离器 Iris Spatula	52404A	短角形、头长11mm、全长 120mm Short angled, 11mm in tip length
7	烧灼止血器 Cauter	52700	球针型、圆形Φ5.5mm、全 长119mm Ball-shape tip, Φ5.5mm, 119mm in total length
8	维纳斯剪 Vannas scissors	54147B	弯形、刃长9mm、全长84mm Curved, 9mm in blade length, 84mm in total length
9	眼用显微手术剪 Micro Ophthalmic Scissors	54022	直圆头、全长85mm Curved round pointed tips, 85mm in total length
10	小梁剪 Micro-corneal Scissors	54108A/B	弯形、刃长12mm、全长120mm Curved, 12mm in blade length, 120mm in total length

序号	产品名称	代号	形式规格
11	缝线结扎镊 Suturing Forceps	53325A	直平台、唇头钩0.15mm、 全长105mm Suturing Forceps Straight, 1x2 teeth of 0.15mm in total length
12	小梁咬切器 Micro trabecula puncher	56105A	Φ1mm、锁式半圆 Φ1mm, semicircle, forceps-type
13	小梁切开刀 Trabeculotomy probe	56092	右式、Φ0.28mm、全长80mm Right Φ0.28mm, 80mm in total length
14	小梁切开刀 Trabeculotomy probe	56090	左式、Φ0.28mm、全长80mm Left Φ0.28mm, 80mm in total length
15	眼用规 Castroviejo Caliper	56211A	可调式测量范围0-20mm Adjustable measuring from 0-20mm at 1mm increment

SYX16 白内障小切口手术器械包

SYX16 Cataract Small Incision Instruments Set



序号	产品名称	代号	形式规格
1-2	隧道刀 Discission Knife	52274A	尖头,刃长9mm,全长128mm Pointed, 9mm in blade length, 128mm in total length
		52275A	圆头,刃长9mm,全长128mm Round, 9mm in blade length, 128mm in total length
3-4	碎核托板 Nucleus Spatula	52334	左式,头部宽2.8mm,全长125mm Left, 2.8mm in head width, 125mm in total length
		52335	右式,头部宽2.8mm,全长125mm Right, 2.8mm in head width, 125mm in total length
5	碎核刀 Cataract Knife	52276D	角形,刃长10.5mm,全长120mm 10.5mm in blade length, 120mm in total length
6	移核镊 Capsul Extractive Forceps	53328A	角形唇头齿距0.5mm,全长105mm Angled, 0.5mm in serration length, 105mm in total length
7	撕囊镊 Capsulorhexis Forceps	53317A	弧形圆头宽0.3mm,全长105mm Round tips, 0.3mm in tip width, 105mm in total length
8	注水手柄 Infusion Handpiece	56203	双接头,全长100mm 100mm in total length
9	白内障注吸器 Irrigation/aspiration Cannula	56191A	双弯管,φ0.7mm φ0.7mm Double curved tubes.
10	晶体线环 Lens Loop	56136A	注水式有齿3.5×11mm 3.5mm(width)×11mm(length) slightly

序号	产品名称	代号	形式规格
11-12	十二点位置注吸器 12 O'clock Cannula	56195	左式, φ0.7mm Left φ0.7mm
		56196	右式, φ0.7mm Right φ0.7mm
13	人工晶体植入镊 IOL Implantation Forceps	53580A	弯形, R12, 全长112mm Curved, R12, 112mm in total length
14-16	膜状内障剪 Cataract Scissors Membrane	54268A/B	直形,刃长16mm,全长120mm Straight, 16mm in blade length, 120mm in total length
		54269A/B	弯形,刃长16mm,全长120mm Curved, 16mm in blade length, 120mm in total length
15	显微持针钳 Micro-needle Holder	55152A/B	弯形,全长110mm Curved, 110mm in total length

SYX17 超声乳化显微手术器械包

SYX17 Microsurgical Instruments Set for Phaco



序号	产品名称	代号	形式规格
1	小梁剪 Micro-corneal Scissors	54108A	弯形,刃长12mm,全长120mm Curved, 12mm in blade length, 120mm in total length
2	开睑器 Eye Speculum	51286A	钢丝式,扩张杆长14mm,全长50mm Wire speculum, stainless 14mm blade in length, 50mm in total length
3	晶体植入切开刀 Discission Knife for IOL Implantation	52279A	角形,刀宽4mm,40°,全长113mm 4mm in blade width, 40°, 113mm in total length
4	角膜膜切开刀 Corneal Membrane Knife	52286A	直形,斜角10°,全长111mm Straight, 10°, 111mm in total length
5	巩膜隧道刀 Scleral Discission Knife	52277A	月形,刀宽2mm,40°,全长113mm 2mm in width of knife, 40°, 113mm in total length
6	裂隙穿刺刀 Slit Knife	52278A	角形,刀宽3mm,40°,全长113mm 3mm in width of knife, 40°, 113mm in total length
7	人工晶体固定钩 IOL Fixation Hook	51474A	角形,φ0.2mm,全长110mm Angled φ0.2mm, 110mm in total length
8	核转动器 Lens Manipulator	51470A	叉形长头宽0.7mm,全长115mm Y-shaped tip, 0.7mm, 115mm in total length
9	核抛光器 Capsule Polisher	51472A	角形,头宽0.6×0.8mm,全长113mm 0.6×0.8mm elliptical tip, 113mm in total length
10	晶体切核器 Lens Nucleus Manipulator	51473A	角形头宽0.25mm,全长112mm 0.25mm in tip width, 112mm in total length

序号	产品名称	代号	形式规格
11	虹膜复位分离器 Iris Spatula	52417A	长角形头宽0.2mm,全长112mm φ0.2mm in tip width, 112mm in total length
12	移核器 Lens Nucleus Hook	51471A	弯角形头宽0.2mm,全长115mm 0.2mm in head width, 115mm in total length
13	白内障注吸器 Aspiration/Irusion of Cataract	56188A	弯三接头,φ0.6mm Curved tip, triple joint, φ0.6mm
14	撕囊镊 Capsulorhexis Forceps	53317A	弧形圆头宽0.3mm,全长105mm Round tips, 0.3mm in tip width, 105mm in total length
15	系结镊 Tying Forceps	53320A	直平台长4.5mm,全长105mm Straight, platform 4.5mm in length, 105mm in total length
16	缝线结扎镊 Suturing Forceps	53325A	直平台长4.5mm,0.15唇头钩,全长105mm Straight, platform 4.5mm in length, 1×2 teeth of 0.15mm, 105mm in total length
17	注水器 Infusion Cannula	56197A	角形,头宽0.5×1.1mm Angled, 0.5×1.1mm elliptical tip

SZY-CBM19眼科手术器械包

SZY-CBM19 Ophthalmic Operation Instrument Set




序号	产品名称	代号	形式规格
1	烧灼止血器 Cauter	52700	球针型、圆形 $\phi 5.5\text{mm}$ 、全长119mm Ball-shape tip, $\phi 5.5\text{mm}$, 119mm in total length
2	维纳斯剪 Vannas scissors	54136B	直形、刃长6.5mm、全长84mm Straight, 6.5mm in blade length, 84mm in total length
3	眼用手术剪 total length	54014	弯圆头、全长100mm Curved blunt tips, 100mm in total length
4	虹膜剪 Iris Scissors	54060	双尖、全长90mm Two sharp tips, 90mm in total length
5	眼用镊 Ophthalmic Forceps	53076	直、1 $\times$ 2唇头钩、全长100mm Straight, 1 $\times$ 2 teeth, 100mm in total length
6	角膜剪 Corneal scissors	54109A/B	弯形、全长120mm Curved, platform 4.5mm in length, 105mm in total length
7	显微持针钳 Micro-needle Holder	55152A/B	弯形、全长110mm Curved, 110mm in total length
8	斜视钩 Strabismus Hook	52454	弯圆头、头宽12mm、全长135mm Blunt, 12mm in tip width, 135mm in total length
9	刀片夹持器 Blade Holder	56202B	弹簧式、全长120mm With spring, 120mm in total length
10	虹膜复位器 Iris Spatula	52391	弯形、头宽1.7mm、全长130mm Curved, 1.7mm in tip width, 130mm in total length

序号	产品名称	代号	形式规格
11	晶体线环 Lens Loop	52500	弯形、头部4.5mm $\times$ 15mm、全长140mm Curved, tip 4.5mm $\times$ 15mm, 140mm in total length
12	开睑器 Eye speculum	51286A	钢丝式、扩张杆长14mm、全长50mm Wire speculum, stainless 14mm blade in length, 50mm in total length
13	无损伤镊 Corneal Scleral Forceps	53286	直半圆、全长70mm Straight, semicircle tips, 70mm in total length
14	虹膜镊 Iris Forceps	53004	曲1 $\times$ 2唇头钩、全长70mm Bend 1 $\times$ 2 teeth, 70mm in total length
15	白内障注吸器 Coaxial I/A Cannula	56191A (含56199)	弯双管、 $\phi 0.7\text{mm}$ Double curved tubes $\phi 0.7\text{mm}$
16	系结镊 Tying Forceps	53322A	弯平台长4.5mm、全长105mm Curved, platform 4.5mm in length, 105mm in total length
17	系结镊 Tying Forceps	53320A	直平台长4.5mm、全长105mm Straight, platform 4.5mm in length, 105mm in total length
18	晶体囊镊 Capsulotomy Forceps	53274	弯长方匙、匙径1.4mm $\times$ 2.4mm、全长70mm Oblong spoon tips, curret 1.4mm $\times$ 2.4mm, 70mm in total length
19	固定镊 Fixation Forceps	53340	直形、齿距0.7mm、全长110mm Straight, 0.7mm in serrated tips distance, 110mm in total length

SYX20 人工晶体白内障显微手术器械包

SYX20 Micro Instrument Set for Cataract and Intraocular Lens Implantation Surgery



序号	产品名称	代号	形式规格
1	白内障注吸器 Coaxial I/A Cannula	56190	直三棱式、 $\phi 0.6\text{mm}$ Straight tip, triple joint, 0.6mm
2-4	冲洗针 Aspiration Needle	56146	角形、 $\phi 0.4\text{mm}$ 、15°、全长27cm Angled $\phi 0.4\text{mm}$ tip, 15°, 27mm in total length
		56144	直形、 $\phi 0.4\text{mm}$ 、全长27mm Straight $\phi 0.4\text{mm}$ tip, 27mm in total length
		56148	直形、 $\phi 0.7\text{mm}$ 、15°、全长27mm Straight $\phi 0.7\text{mm}$ tip, 15°, 27mm in total length
5	白内障注吸器 Coaxial I/A Cannula	56191A	弯双管、 $\phi 0.7\text{mm}$ Double curved tubes, $\phi 0.7\text{mm}$
6	开睑器 Eye Speculum	51288A	钢丝式、扩张杆长15mm、全长50mm Wire speculum, stainless 15mm, blade in length, 50mm in total length
7	刀片夹持器 Blade Holder	56202B	弹簧式、全长120mm With spring, 120mm in total length
8	白内障匙 Lens Spoon	52574	匙径8 $\times$ 2.5mm、全长130mm 8 $\times$ 2.5mm in spoon length, 130mm in total length
9	晶体线环 Lens Loop	52488A	弯形、头部4 $\times$ 6mm、全长120mm Tip 4mm $\times$ width $\times$ 6mm (length) curved loop, 120mm in total length
10	虹膜复位分离器 Iris Spatula	52404A	短角形、头长11mm、全长120mm Curved, 11mm in tip length, 120mm in total length

序号	产品名称	代号	形式规格
11-12	人工晶体定位钩 IOL Fixation Hook	51474A	角形、 $\phi 0.2\text{mm}$ 、全长110mm Angled $\phi 0.2\text{mm}$, 110mm in total length
		51478A	弯形、 $\phi 0.2\text{mm}$ 、15°、全长110mm 0.2mm cylindrical tip, angled 15°, 110mm in total length
13-14	维纳斯剪 Vannas Scissors	54138B	直形、刃长6.5mm、全长90mm Straight, 6.5mm in blade length, 90mm in total length
		54140B	弯形、刃长6.5mm、全长84mm Curved, 6.5mm in blade length, 84mm in total length
15	显微持针钳 Micro-needle Holder	55152B	弯形、全长110mm Curved, 110mm in total length
16	系结镊 Tying Forceps	53320A	直平台长4.5mm、全长105mm Straight, platform 4.5mm in length, 105mm in total length
17	角膜剪 Corneal Scissors	54109B	弯形、刃长14mm、全长120mm Slight curved, 14mm in blade length, 120mm in total length
18	缝线结扎镊 Suturing Forceps	53325A	直平台长4.5mm、0.15唇头钩、全长105mm Straight, platform 4.5mm in length, 1 $\times$ 2 teeth of 0.15mm, 105mm in total length
19	膜状白内障剪 Cataract Scissors Membrane	54264B	直形、刃长10mm、全长117mm Straight tips, 10mm in blade length, 117mm in total length
20	系结镊 Tying Forceps	53324A	角平台长4.5mm、全长105mm Short angled, platform 4.5mm in length, 114mm in total length

SZY-CBM21眼科手术器械包

SZY-CBM21 Ophthalmic Operation Instrument Set



序号	产品名称	代号	形式规格
1	固定镊 Fixation Forceps	53340	直形,齿距0.7mm,全长110mm Straight, 0.7mm in serrated tips distance, 110mm in total length
2-3	晶体囊镊 Capsulotomy Forceps	53274	弯长方匙,匙径1.2mm×3mm,全长70mm Oblong spoon tips, curette 1.2mm×3mm, 70mm in total length
4	系结镊 Tying Forceps	53323A	角平台,长4.5mm,全长105mm Long angled platform 4.5mm in length, 105mm in total length
5	系结镊 Tying Forceps	53320A	直平台,长4.5mm,全长105mm Straight platform 4.5mm in length, 105mm in total length
6-7	虹膜镊 Iris Forceps	53042	弯1×2唇头钩,全长60mm Curved 1×2 teeth, 60mm in total length
		53004	曲1×2唇头钩,全长70mm Bend 1×2 teeth, 70mm in total length
8	无损伤镊 Corneal Scleral Forceps	53286	直半圆,全长70mm Straight, Semicircle, 70mm in total length
9-10	开睑器 Eye Speculum	51286A	铰链式,扩张杆长14mm,全长50mm Open blades, stainless 14mm blade in length, 50mm in total length
		51262	铰链式,扩张杆长14mm,全长75mm Open blades, stainless 14mm blade in length, 75mm in total length
11	晶体线环 Lens Loop	52500	弯形,头部4.5mm×15mm,全长140mm Curved tip 4.5mm (width)×15mm (length), 140mm in total length

序号	产品名称	代号	形式规格
12	斜视钩 Strabismus Hook	52454	弯圆,头宽12mm,全长13.5mm Blunt, curved tip, 12mm in head width, 13.5mm in total length
13	刀片夹持器 Blade Holder	56202B	弹簧式,全长120mm Spring type, 120mm in total length
14	虹膜复位器 Iris Spatula	52391	弯形,头宽1.7mm,全长130mm Curved, 1.7mm in tip width, 130mm in total length
15-16	显微持针钳 Micro-needle Holder	55152A/B	弯形,全长110mm Curved, 110mm in total length
		55115	直小方柄,全长130mm Straight, 130mm in total length
17	眼用镊 Ophthalmic Forceps	53076	直1×2唇头钩,全长100mm Straight, 1×2 teeth, 100mm in total length
18	虹膜剪 Iris Scissors	54060	双尖,全长90mm Double pointed tips, 90mm in total length
19	眼用手术剪 Ophthalmic Scissors	54006	直圆头,全长100mm Straight blunt tips, 100mm in total length
20	维纳斯剪 Vannas Scissors	54136B	直形,刃长6.5mm,全长84mm Straight, 6.5mm in blade length, 84mm in total length
21	角膜剪 Corneal Scissors	54102	弯形,R5.5mm Curved, R5.5mm

SYX24斜视手术器械包

SYX24 Strabismus Surgical Instruments Set



序号	产品名称	代号	形式规格
1	斜视镊 Strabismus Forcep	53801A	右弯平台, 0.25唇头钩 Right curved platform, 1×2 teeth of 0.25mm
2	斜视镊 Strabismus Forcep	53800A	左弯平台, 0.25唇头钩 Left curved platform, 1×2 teeth of 0.25mm
3	斜视镊 Strabismus Forcep	53802A	直式平台, 0.25唇头钩 Straight platform, 1×2 teeth of 0.25mm
4	系结镊 Tying Forcep	53320A	直平台,长4.5mm, 全长105mm Straight platform 4.5mm in length, 105mm in total length
5	斜视镊 Strabismus Forcep	53583	直式平台, 0.35唇头钩 Straight platform, 1×2 teeth of 0.35mm
6	显微持针钳 Micro Needle Holder	55159T	弯锁止钩, 全长140mm Curved with lock, 140mm in total length
7	结膜剪 Conjunctival Scissor	54095A	弯形,刃长17mm, 全长115mm Curved, 17mm in blade length, 115mm in total length
8	小梁剪 Micro Corneal Scissor	54108B	弯尖头,刃长12mm, 全长120mm Curved round tip, 12mm in blade length, 120mm in total length
9	开睑器 Eye Speculum	51280	铰链式,板宽6mm,扩张杆长14mm, 全长64mm Solid Blades, 6mm in board width, 14mm blades in length, 64mm in total length
10	开睑器 Eye Speculum	51280A	铰链式,板宽8mm,扩张杆长14mm, 全长64mm Solid Blades, 8mm in board width, 14mm blades in length, 64mm in total length
11	斜视钩 Strabismus Hook	52440	弯圆头,头宽1.6mm Curved round tip, 1.6mm in tip width
12	斜视钩 Strabismus Hook	52439A	弯圆头,头宽1.8mm Curved round tip, 1.8mm in tip width

序号	产品名称	代号	形式规格
13	斜视钩 Strabismus Hook	52457A	立体钩,椭圆头0.7mm 3dimensional hook, elliptical head 0.7mm
14	斜视钩 Strabismus Hook	52452A	弯圆头,头宽1.3mm Curved round tip, 1.3mm in tip width
15	斜视钩 Strabismus Hook	52451A	弯圆头,头宽1.2mm Curved round tip, 1.2mm in tip width
16	斜视钩 Strabismus Hook	52459A	F钩头,头宽1.5mm F head, 1.5mm in tip width
17	斜视钩 Strabismus Hook	52456A	钝钩头,头宽2.5mm Dull tip, 2.5mm in tip width
18	斜视钩 Strabismus Hook	52455A	弯圆头,Φ0.6mm Curved round tip, Φ0.6mm
19	斜视钩 Strabismus Hook	52460A	双钩头,Φ0.8mm Double head, Φ0.8mm
20	眼睑拉钩 Eyelid Hook	51452A	头宽10mm 10mm in tip width
21	斜视圆挑板 Strabismus Hook	52462	头宽7.3mm 7.3mm in tip width
22	斜视圆挑板 Strabismus Hook	52462A	头宽9.2mm 9.2mm in tip width
23	斜视圆挑板 Strabismus Hook	52463	头宽6mm 6mm in tip width
24	斜视钩 Strabismus Hook	52458A	刺腹头,头宽1.75mm Head with dial, 1.75mm in tip width